



POST-OPERATIVE AND CARE INSTRUCTIONS FOR VISICLEAR/VALPLAST/FLEXIBLE REMOVABLE PARTIAL DENTURE

- Clean denture with water or you can use denture soaking tablets (i.e. Efferdent, Polygrip, generic brand). Soak the partial for 5-10 minutes in warm water then rinse.
- When not wearing your denture, store in a wet container.
- Make sure to remove your denture prior to brushing your teeth. Toothbrushes and tooth paste contain materials that may scratch and damage your denture.
- Please keep your denture away from pets, especially dogs. A case is being provided for you, and we strongly suggest when not wearing it, store it in the case.
- Do not bend or twist your Flexible partial denture.
- After having time to wear our appliance, you may notice ulcerations/sore spots. If this occurs, please call our office for an adjustment (at no charge).
- It can take up to 3 months for your appliance to become a natural part of your body, be persistent! During the adjustment period start with small bites and slow eating.

*Thank you,
Dr. Gupta, and Staff*

2800 Niles Road, Saint Joseph, MI 49085

Tel: 269-429-2511 | Fax: 269-429-5130 | info@smilesonniles.com | www.smilesonniles.com